

In-Use Storage and Handling of Large Cylinders – Streamlined

Time saving, easy to load, easy to move carts minimize in-room footprint for large cylinders

The new Praxair *Medipure*[™]-*LC* Single and Dual Cylinder Carts makes having cylinders ready to go at the point of care easier than ever. Available in either a two cylinder or a single cylinder format, these advanced carts allowing right storage of cylinders and reducing cylinder footprint in the ER or ICU areas. *Medipure-LC Series* carts improve cylinder handling and make in-room storage more functional and aesthetically pleasing.

Features	Benefits
Single and Dual Cylinder Design	Your choice of formats to maximize cylinder storage volumes and footprint to match your space constraints
Vertical Storage Cylinder Position	Replaces traditional tube-metal cylinder carts which require storage of cylinders at an angle for stability
Accessory Baskets	Store hoses, vinyls and other items in a neat ready- to-use basket Eliminates the need to 'hang' items from the cylinder and connections
Sturdy Construction	Robust design facilitates storage of 'K/H,' 'T,' or '6K' cylinders
Smooth Easy Wipe-Down Surfaces	Smooth surfaces make it easier to facilitate care and cleaning to help maintain cleanliness in patient care areas
Easy Lift Mechanism	Crank-type design reduces the effort required get cylinders into a position that's ready for transport
Locking Wheels	The cart is effectively held into position
IV Pole Mount	Mount gas conditioning and other devices required for patient care and maximize space utilization
Optional Cylinder mount- ing bracket for portable gas cylinders*	Portable medical gas cylinders in 'D' and 'E' size can be mounted when required for transport



^{*}Available for Single Cylinder Cart only









Single Cylinder Design

OR – High Pressure N₂